

Beltone Direct Remote Control 2

Beltone continues to lead the hearing instrument industry in wireless technology innovation. 2.4 GHz technology allows clients to receive clear audio signals from a variety of sources directly to their hearing aids. No neckloops required.

Now, the Beltone Direct Remote Control 2 gives Beltone hearing aid wearers even greater control over their very own Personal Wireless Network.



What's New with the Direct Remote Control 2?

- **Adjustable microphone and audio streaming balance** to make sure that the audio from a desired source and sounds in the surrounding environment are both at comfortable volume levels.
- **Mute button** to hear only the audio source of interest when surrounding noise becomes too loud.
- **New antenna design** for stronger connection to hearing instruments and uninterrupted use.
- **Color display** with new, easy-to-read icons.
- **Rechargeable battery** for convenient use without worrying about changing batteries.



Beltone **2.4 GHz Personal Wireless Network**

Worldwide Headquarters

Beltone A/S
Lautrupbjerg 7
DK-2750 Ballerup
Denmark
Tel.: + 45 45 75 11 11
Fax: + 45 45 75 11 19
beltone-hearing.com

