







# SolusPro Cable Guide

Circuit Model	 <b>Grey Cable:</b> Cable # SA04305 Red Cable # SA04304 Blue Black Boot # 18041400 Taupe Boot # 17598200	 <b>Grey Cable:</b> Cable # SA04305 Red Cable # SA04304 Blue Flex Strip # 19533900	 <b>Grey Cable:</b> Kit Part # 19215000	<b>IIC</b>  <b>Battery Pill Cables:</b> Part # 18523700 Part # 18523800	<b>Prime+/Prime</b>  <b>Grey Cable:</b> Cable # SA04305 Red Cable # SA04304 Blue Battery Pill # 17167800	<b>Obsolete</b>  <b>Black Cable:</b> Cable # SA04302 Flex Strip # 1730901
<b>Receiver in the Ear (RIE)</b>						
Legend, Ally, Promise, True & Origin 64	● Black					
Legend 63		●				
First, Promise, True 63		●				●
Promise and True 62		●				●
<b>Behind the Ear (BTE)</b>						
Legend & Ally 86, 76	● Black					
Ally 66	● Black					
Promise & True 78, 66	● Black					
Promise 85	● Black					
Origin 85, 75, 65 Bold 95	● Taupe					
Boost 95			●			
<b>Custom Products</b>						
All Legend & Ally		●				
All Promise, True & Origin		●				●
micro-Invisa/IIC				●		
Prime+ and Prime					●	