

Beltone App Compatibility Chart

	Apple (MFi direct streaming)	Android ** (Phone Link 2 needed for audio streaming)
Amaze, Trust, Boost Max <u>HearMax</u> (iOS 10 or higher, OS 5.0 Lollipop or higher) (Watch OS 3)	iPhone: X-5 th generation iPad: Pro, Air2, Air, mini 4-2, iPad 6-4 th generation iPod touch 6 th generation (iOS 8 or later) Apple Watch: 1-2 series	Samsung Galaxy: S9+-S4, Note 5-4 Sony Xperia Z3 Google Pixel XL, Doro 8030 HTC: 10, One M8, One M9 LG: G3, G5, V20 Huawei: Mate M8, P9
Legend, Boost Plus, First, Boost <u>HearPlus</u> (iOS 8 or higher, OS 5.0 Lollipop or higher) (Watch OS 2 or higher w/ iOS 9 or higher)	iPhone: X-5 th generation iPad: Pro, Air2, Air, mini 4-2, iPad 6-4 th generation iPod touch 5-6 th generation Apple Watch: 1-2 series	Samsung Galaxy: S9+-S4, Note 5-4 & 8
Promise, Ally 3 (needs Phone Link2 for app to work) Smart Remote	Any	Any Android phone
Any product Tinnitus Calmer	Any	Any Android phone

Apple, the Apple logo, iPhone, iPad Pro, iPad Air, iPad mini, iPad and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries.

**We recommend that your Android mobile device runs Android OS version 5.0 - Lollipop or later. Your hearing aids may need to be updated to the latest firmware to work with the app. This update can be done by your hearing care professional.

For a full list of compatible phones and the latest phone additions please go to:

<http://www.beltone-hearing.com/en/support/compatibility>