



Product Feature Comparison

-  Ultimate
-  Advanced
-  Basic

	Beltone Imagine 17	Beltone Imagine 9	Beltone Imagine 6
Sound Quality & Hearing in Noise			
Directional Options for Better Hearing in Noise	 CrossLink Directionality 3	 Crosslink Directionality 2 w/Personal Sound ID	 Crosslink Directionality
M&RIE for naturally shaped sound			
Improved spacial awareness			
Binaural Beamformer for improved speech recognition in noise		 Only when using M&RIE	
Band Split Directionality for Improved Sound Quality	 With adjustable mixing point	 With adjustable mixing point	
Ultra Focus for pointed directionality in noisy environments			
Comfort in Noise			
Per Environment Sound Cleaner Option for Automatic Noise Reduction			
Sound Cleaner for Noise Reduction	 Mild, Moderate, Considerable, Strong	 Mild, Moderate, Considerable	 Mild and Moderate
Impulse Noise Reduction to Increase Comfort for Loud Noises			 Mild and off
Wind Noise Reduction	 Mild, Moderate, Strong	 Mild, Moderate,	 On/Off
Personalization			
Smart Gain for Automatic Volume Adjustments			
WhistleStop Feedback Eraser Setting for Maximum Feedback Control			
Feedback Eraser	 Mild, Moderate, Strong, Music Mode	 Mild, Moderate, Strong, Music Mode	 Mild, Moderate, Music Mode
HearMax App Sound Enhancer Controls for Maximum Personalization			